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[www.deep-company.com](http://www.deep-company.com)

[contact@deep-company.com](mailto:contact@deep-company.com)

**IMPORTANT : Please read these instructions carefully before use**

### LIABILITY

The user of DEEP Company ancillary components has the duty to determine whether or not any product is suitable for the particular patient and circumstances. DEEP Company disclaims any liability, express or implied, and shall have no responsibility for any, direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of DEEP Company products. The clinician is also obliged to study the latest developments in regard to the DEEP Company products and their applications regularly. In cases of doubt, the clinician has to contact DEEP Company. Since the processing and surgical application of this product are under the control of the clinician, they are his/her responsibility. DEEP Company assumes no liability whatsoever for damages arising thereof.

### INDICATIONS/CONTRAINDICATIONS

Each DEEP Company ancillary has been designed to meet specific indications: Refer to the IFU of the implant or the associated prosthetic system, and the specific IFU supplied with the ancillary if necessary. DEEP Company disclaims any liability for non-compliance with indications and contraindications of these devices.

### INSTRUCTIONS FOR CLINICIAN

Each implant system has specific design characteristics for abutments and prosthetic components. Combining components that are not dimensioned or configured for correct mating can lead to mechanical failure of components, tissue damage or unsatisfactory aesthetic results.

We strongly recommend that clinicians, new as well as experimented implant users, always go through special training before undertaking a new treatment method. Always work with an experienced colleague the first time you employ a new treatment method.

### PROCEDURAL PRECAUTIONS

The responsibility of the practitioner :

- To prevent the risk of inhalation and ingestion by the patient.
- To make sure to have the necessary ancillary for handling implant or prosthetic system.
- To check the screw torques, ancillaries functioning and retention force before use.
- To use the shortest drill as possible in accordance with the clinical requirements.
- To externally irrigate the surgical site while drilling.
- To not file contaminated instrument(s) with clean instrument(s) in the same container.

The drilling speed should not exceed 2000 rev / min.

The implant screwing speed should not exceed 25rev / min.

We recommend 20 uses for drills (drilling or sterilization).

Warning : Any incident related to ancillary should be reported to DEEP Company.

### POTENTIAL ADVERSE EFFECTS

- Dislocation, rupture or loss of a component, infections, aesthetic problems, allergies to the different materials contained in prosthetic components.
- Aspirations of a component by the patient, allergy to any component
- Transient or permanent local nervous lesion, local sensitivity loss
- Bone fracture
- Infection

The surgeon is responsible for the complications that can result from an erroneous indication, of a faulty surgical technique, or a lack of sterility. In any case, these complications cannot be attributed to DEEP Company.

### PATIENT INFORMATION:

The clinician is responsible for:

- informing the patient of the potential risks and adverse effects associated with the fitting of such an implant system,
- informing the patient that the effectiveness of the devices depends on his behaviour and his oral hygiene,
- ensuring that the patient reports him all suspicious movements between implant and abutment, and will have to undergo periodic inspections for detecting signs that precede abnormal functional symptoms.

### PACKAGING, STORAGE AND HANDLING

From delivery until use, store devices in a clean room away from products which could have a corrosive action.

In the event of an apparent packaging defect, the product must be returned to DEEP Company.

The non-sterile equipment supplied by DEEP Company has been cleaned of manufacturing residues and then washed, but does not have a sufficient level of decontamination to permit it to be used directly. A decontamination and cleaning treatment is therefore essential prior to any sterilization.

Avoid contact with any objects, liquids or chemicals that could alter the surface of the prosthesis system.

### PREPARATION OF DEVICE FOR AUTOCLAVE STERILIZATION

Composition of detergents and/or disinfectants: we recommend the use of neutral detergent for the cleaning of DEEP Company products. Acid and alkaline detergents should be avoided as well as the use of dishwasher shelves or any household products.

Avoid letting any soil dry onto reusable instruments after they have been used. Every medical device must undergo pre-disinfection (soaking in an aldehyde-free detergent-disinfectant bath). Prolonged contact with blood, pulp, bone, dentine, cement and other residues could damage the devices and compromise the effectiveness of the entire treatment.

### LIMITATIONS ON REPROCESSING OF REUSABLE MEDICAL DEVICES

Repeated cleaning and sterilisation cycles have little effect on medical devices. The end-of-life for a reusable medical device is usually determined by the amount of wear and damage from a use or reprocessing different from the one recommend by DEEP Company.

### PROTOCOL BEFORE USE

<b>Place of use</b>	Do not return soiled medical devices to their storage box. Disassemble any medical devices having multiple parts. Remove any significant residue on medical devices by rinsing them with cold water or wiping them with lint-free tissue.
<b>Containment &amp; transport</b>	It is imperative to reprocess medical devices immediately after use

<b>Manual Cleaning</b>	<p><b>Equipment:</b> Detergent (DENTASEPT ENZYMATIQUE d'Anios - Concentration 1%), non-metallic soft brush, tap water, demineralised water (flushing)</p> <p><b>Note:</b> Be careful not to damage medical devices. Do not bang them together.</p> <ol style="list-style-type: none"> <li>1. Fully submerge the devices in a cleaning bath having a temperature of less than 50°C (Room temperature &lt; 30°C). Any hinged devices should be fully opened and any multicomponent devices should be disassembled. Follow the soaking time recommended by the product manufacturer (15 minutes).</li> <li><b>Note:</b> Unless specified otherwise, instruments and devices must be soaked for at least 15 minutes but avoid over soaking.</li> <li>2. Brush the medical device inside and out, and any engraving, to remove persistent stains. Run the solution through device cavities to remove any residue.</li> <li>3. Rinse thoroughly (Temp &lt; 30°C) by soaking and/or spraying with water to remove detergent and soil; demineralised water should be used whenever possible.</li> </ol>
<b>Manual disinfection</b>	<p><b>Note:</b> Some detergents also act as disinfectants (see manual cleaning step). Soak in the disinfectant solution (Dentasept® Enzymatique from Anios - Concentration 1%) for at least 15 minutes, and then rinse with demineralised water.</p>
<b>Automated cleaning</b>	<p><b>Equipment:</b> Washer-disinfector according NF EN ISO 15883-1&amp;-2 (MIELE G 7825), Detergent (Néodisher® MediClean - concentration 0.5%), rinsing agent (Néodisher® MediKlar - concentration 0.05%), demineralised water.</p> <ol style="list-style-type: none"> <li>1. Arrange the devices in a manner that allows water to flow into their cavities.</li> <li>2. Start the cycle:                     <ul style="list-style-type: none"> <li>- Prewash without product for 1 minute at room temperature</li> <li>- Wash: 40°C water + detergent</li> <li>- Bring temperature up to 50°C and maintain it for 5 minutes</li> <li>- Rinse with cold water for 1 minute at room temperature</li> </ul> </li> <li>3. Unload the instruments and inspect the surface and cavities of each device to make sure all visible soil has been removed. If needed, repeat the cycle or perform manual cleaning.</li> </ol>
<b>Automated disinfection</b>	If using a washer-disinfector, a final rinsing step with thermal disinfection (demineralised water temperature increased to 90°C and held for 5 minutes) can be added. The hold time depends on the size and power of the unit (ISO 15883-1).
<b>Drying</b>	<p><b>Note:</b> This step prevents microbial growth, removes any traces left from the washing and rinsing steps and optimises sterilisation effectiveness.</p> <ul style="list-style-type: none"> <li>- If using the <b>washer-disinfector</b>, do not exceed 120°C (Drying for 20 min at 95°C).</li> <li>- <b>Manual drying</b> can be performed using a clean non-woven absorbent tissue or a clean lint-free cloth, and medical compressed air.</li> </ul>
<b>Packaging</b>	<p><b>Note:</b> No matter which cleaning method is used, medical devices must be quickly placed into packages to avoid recontamination.</p> <p>Use materials suitable for sterilisation; the packaging used to ship the prosthesis systems cannot be sterilised. Make sure the pouch is large enough to hold the medical device, without stretching the edges. Do not place multiple devices in the same pouch.</p>

<b>Sterilisation</b>	Autoclave with high pressure saturated steam at 134°C for at least 18 minutes (prion cycle). If several medical devices are being sterilised at once, make sure the autoclave is not overloaded.
<b>Storage</b>	Isolate the sterile packaged items so they do not come into contact with other devices.

### RETURN:

In the event of device return following an order error, it must be returned unused in its original unopened packaging.

For any device return to after-sales service, please clean and sterilize the device(s) according to the previous described procedures and send them in double sterile pouch.

Product Family	Intended use	Contra-indications	Material / Handling
Drills	Drill into the jaw's bone	<ul style="list-style-type: none"> <li>• Drill into another body part</li> <li>• Use with a turbine</li> </ul>	Stainless Steel Contra-angle connexion
Contra-angle wrenches	Screw and unscrew	Use with a turbine	Stainless Steel Contra-angle connexion
Implant Driver	Assembly compatibility	<ul style="list-style-type: none"> <li>• Use with incompatible implant</li> <li>• Rotation speed or tightening torque inappropriate</li> </ul>	Stainless Steel Each driver is compatible with only one implant range
Drill extension	Facilitate the drill's use between teeth	<ul style="list-style-type: none"> <li>• Insufficient oral aperture</li> </ul>	Stainless Steel Contra-angle connexion
Axis indicator	Control the drilling axis	NA	Grade 5 titanium Manual use

<b>SYMBOLS</b>	 See to instruction for use	<b>LOT</b> Batch number	 0123 EC mark delivered by Notified Body	 Do not re-use
----------------	--	-------------------------	---	---

**DEEP Company thanks you for your support and remains at your disposal for any additional information you may require**

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