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[www.deep-company.com](http://www.deep-company.com)

contact@deep-company.com

**IMPORTANT : Please read these instructions carefully before use**

### LIABILITY

The user of DEEP Company prosthetic components has the duty to determine whether or not any product is suitable for the particular patient and circumstances. DEEP Company disclaims any liability, express or implied, and shall have no responsibility for any, direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of DEEP Company products. The clinician is also obliged to study the latest developments in regard to the DEEP Company products and their applications regularly. In cases of doubt, the clinician has to contact DEEP Company. Since the processing and surgical application of this product are under the control of the clinician, they are his/her responsibility. DEEP Company assumes no liability whatsoever for damages arising thereof.

### DESCRIPTION

Premanufactured endosseous dental prosthetic parts directly connected to the dental implant, intended for use as an aid in the prosthetic rehabilitation.

### INDICATIONS

Prosthetic components are used to make prosthetic rehabilitations that range from replacing one single tooth to an entire arch of bridgework, as well as retentive elements for overdenture applications.

### CONTRA-INDICATIONS

The prosthetic design must accommodate individual patient conditions such as bruxism or unfavorable jaw relationships to reduce the risk of overload or unusual wear, and treatment is contraindicated if adequate accommodation cannot be accomplished.

Temporary plastic abutments are only intended to be used for maximum 90 days in the oral cavity. For titanium abutments, no specific time limit applies.

### INSTRUCTIONS FOR CLINICIAN

Each implant system has specific design characteristics for abutments and prosthetic components. Combining components that are not dimensioned or configured for correct mating can lead to mechanical failure of components, tissue damage or unsatisfactory aesthetic results.

We strongly recommend that clinicians, new as well as experimented implant users, always go through special training before undertaking a new treatment method. Always work with an experienced colleague the first time you employ a new treatment method.

### PROCEDURAL PRECAUTIONS

Close cooperation between surgeon, restorative dentist and dental laboratory technician is essential for a successful implant treatment.

It is especially important to achieve proper stress distribution through passive adaptation and fitting of the bridge to the implant abutments, adjusting occlusion to the opposing jaw; avoiding excessive transverse loading forces, particularly in immediate loading cases. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient.

### PATIENT INFORMATION:

The clinician is responsible for:

- informing the patient of the potential risks and adverse effects associated with the fitting of such an implant system,
- informing the patient that the effectiveness of the devices depends on his behaviour and his oral hygiene,
- ensuring that the patient reports him all suspicious movements between implant and abutment, and will have to undergo periodic inspections for detecting signs that precede abnormal functional symptoms.

### INTERACTION WITH MAGNETIC RESONANCE IMAGING AND OTHER ELECTROMAGNETIC SYSTEM

The raised temperature after interaction between MRI and implants does not pose any known hazard to the patient. But the quality of the MR image may be compromised, if the area of interest is in the same area of the device: an artefact may appear on the picture.

### POTENTIAL ADVERSE EFFECTS

Dislocation, rupture or loss of a component, infections, aesthetic problems, allergies to the different materials contained in prosthetic components.

### PACKAGING, STORAGE AND HANDLING

From delivery until use, store devices in a clean room away from products which could have a corrosive action.

In the event of an apparent packaging defect, the product must be returned to DEEP Company.

The non-sterile equipment supplied by DEEP Company has been cleaned of manufacturing residues and then washed, but does not have a sufficient level of decontamination to permit it to be used directly. A decontamination and cleaning treatment is therefore essential prior to any sterilization.

Avoid contact with any objects, liquids or chemicals that could alter the surface of the prosthesis system.

### PREPARATION OF DEVICE FOR AUTOCLAVE STERILIZATION

Composition of detergents and/or disinfectants: we recommend the use of neutral detergent for the cleaning of DEEP Company products. Acid and alkaline detergents should be avoided as well as the use of dishwasher shelves or any household products.

Avoid letting any soil dry onto reusable instruments after they have been used. Every medical device must undergo preinfection (soaking in an aldehyde-free detergent-disinfectant bath). Prolonged contact with blood, pulp, bone, dentine, cement and other residues could damage the devices and compromise the effectiveness of the entire treatment.

### PROTOCOL BEFORE USE

<b>Manual Cleaning</b>	<p><b>Equipment:</b> Detergent (DENTASEPT ENZYMATIQUE d'Anios - Concentration 1%), non-metallic soft brush, tap water, demineralised water (flushing)</p> <p><b>Note:</b> Be careful not to damage medical devices. Do not bang them together.</p> <ol style="list-style-type: none"> <li>1. Fully submerge the devices in a cleaning bath having a temperature of less than 50°C (Room temperature &lt; 30°C). Any hinged devices should be fully opened and any multicomponent devices should be disassembled. Follow the soaking time recommended by the product manufacturer (15 minutes).</li> <li><b>Note:</b> Unless specified otherwise, instruments and devices must be soaked for at least 15 minutes but avoid over soaking.</li> <li>2. Brush the medical device inside and out, and any engraving, to remove persistent stains. Run the solution through device cavities to remove any residue.</li> <li>3. Rinse thoroughly (Temp &lt; 30°C) by soaking and/or spraying with water to remove detergent and soil; demineralised water should be used whenever possible.</li> </ol>
<b>Manual disinfection</b>	<p><b>Note:</b> Some detergents also act as disinfectants (see manual cleaning step). Soak in the disinfectant solution (Dentasept® Enzymatique from Anios - Concentration 1%) for at least 15 minutes, and then rinse with demineralised water.</p>
<b>Automated cleaning</b>	<p><b>Equipment:</b> Washer-disinfector according NF EN ISO 15883-1&amp;-2 (MIELE G 7825), Detergent (Néodisher® MediClean - concentration 0.5%), rinsing agent (Néodisher® MediKlar - concentration 0.05%), demineralised water.</p> <ol style="list-style-type: none"> <li>1. Arrange the devices in a manner that allows water to flow into their cavities.</li> <li>2. Start the cycle:                     <ul style="list-style-type: none"> <li>- Prewash without product for 1 minute at room temperature</li> <li>- Wash: 40°C water + detergent</li> <li>- Bring temperature up to 50°C and maintain it for 5 minutes</li> <li>- Rinse with cold water for 1 minute at room temperature</li> </ul> </li> <li>3. Unload the instruments and inspect the surface and cavities of each device to make sure all visible soil has been removed. If needed, repeat the cycle or perform manual cleaning.</li> </ol>
<b>Automated disinfection</b>	<p>If using a washer-disinfector, a final rinsing step with thermal disinfection (demineralised water temperature increased to 90°C and held for 5 minutes) can be added. The hold time depends on the size and power of the unit (ISO 15883-1).</p>
<b>Drying</b>	<p><b>Note:</b> This step prevents microbial growth, removes any traces left from the washing and rinsing steps and optimises sterilisation effectiveness.</p> <ul style="list-style-type: none"> <li>- If using the <b>washer-disinfector</b>, do not exceed 120°C (Drying for 20 min at 95°C).</li> <li>- <b>Manual drying</b> can be performed using a clean non-woven absorbent tissue or a clean lint-free cloth, and medical compressed air.</li> </ul>
<b>Packaging</b>	<p><b>Note:</b> No matter which cleaning method is used, medical devices must be quickly placed into packages to avoid recontamination.</p> <p>Use materials suitable for sterilisation; the packaging used to ship the prosthesis systems cannot be sterilised. Make sure the pouch is large enough to hold the medical device, without stretching the edges. Do not place multiple devices in the same pouch.</p>

<b>Sterilisation</b>	Autoclave with high pressure saturated steam at 134°C for at least 18 minutes (prion cycle). If several medical devices are being sterilised at once, make sure the autoclave is not overloaded.
<b>Storage</b>	Isolate the sterile packaged items so they do not come into contact with other devices.

### EQUIVALENCE CHART:

Our devices ranges are compatible with different implant systems mentioned in the following equivalence chart.

The reference prefix refers to the name of the original range

Prefix	Ranges	Equivalence
BS	LIKE-BS	Brånemark System® from NOBELBIOCARE®
		Universal+® from EUROTEKNICA®
R	LIKE-R	NobelReplace® from NOBELBIOCARE®
A ou CC	LIKE-CC	NobelActive® and ReplaceCC® from NOBELBIOCARE®
S	LIKE-EV	EVL®S and EVL®Konik from GLOBAL D®
ST	LIKE-ST	SynOcta® from STRAUMANN®
I	LIKE-I	InKone® from GLOBAL D®
Z	LIKE-Z	Screw-Vent® from ZIMMER

® Copyrighted Trademark

DEEP Company remains at your disposal for any question relating to the compatibility of their devices with a given implant system.

### RETURN:

In the event of return device, it must be returned unused in its original unopened packaging.

### ESTABLISHMENT-REMOVAL:

Abutment tightening and removal must be performed using appropriate instruments. Screwdriver tips must especially match to the head of the compatible system screw.

### Screw tightening torque:

	LIKE-I	LIKE-EV	LIKE-BS/R/CC et LIKE-ST	LIKE-Z
Direct implant monobloc abutments	20	[30 ;35]	35	30
Direct implant screw-retained abutments	15	[30 ;35]	35	30
Cover and healing screws	Manuel	Manuel	Manuel	Manuel
Indirect screw-retained abutments	15	[30 ;35]	15	20

DEEP Company thanks you for your support and remains at your disposal for any additional information you may require

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