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DEEP Company
ZA du bois Saint-Pierre
38280 JANNEYRIAS - FRANCE

Tel : +33 (0) 4 28 01 26 03
www.deep-company.com
contact@deep-company.com

IMPORTANT : Please read these instructions carefully before use

LIABILITY

The user of DEEP Company ancillary components has the duty to determine whether or not any product is suitable for the particular patient and circumstances. DEEP Company disclaims any liability, express or implied, and shall have no responsibility for any, direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of DEEP Company products. The clinician is also obliged to study the latest developments in regard to the DEEP Company products and their applications regularly. In cases of doubt, the clinician has to contact DEEP Company. Since the processing and surgical application of this product are under the control of the clinician, they are his/her responsibility. DEEP Company assumes no liability whatsoever for damages arising thereof.

DESCRIPTION AND INTENDED USE

LIKE Ancillaries are medical devices allowing the placement of LIKE Implants and prosthetics parts in the context of dental rehabilitations. The intended use of each instrument are given below :

Family product	Intended use
Drill	Drill into alveolar bone
Implant driver (manual and contra angle)	Take the implants in their packaging and place them in the alveolar bone
Drill extension	Facilitate the drill's use between the teeth
Axis indicator	Control the axis of the implant and the depth of preparation
Screwing key / Screwing	Screw/ Unscrew the prosthetic components with the implant
Torque wrench	Control the tightening torque of the screws

Materials :

Family product	Materials
Drill	Stainless steel
Implant driver (manual and contra angle)	
Drill extension	
Screwing key / Screwing	
Torque wrench	
Axis indicator	Titane biocompatible de grade 5

INDICATIONS

LIKE instruments are used during dental surgery for dental rehabilitation in mandible and maxilla. These restorations range from a single replacement to the use of elements for overdenture applications as well as the full arch of bridgework.

INTENDED PATIENT POPULATION

LIKE Implant System is intended for people with the absence of at least one tooth who does not have any contraindications identified in the instructions for use.

For children, the placement of implants is not recommended until the end of jaw growth has been duly established. Preoperative gum and bone deficiencies can compromise the aesthetic result or result in unfavourable implant positioning.

CONTRAINDICATIONS

The contraindications associated with the placement of implants or prostheses by LIKE instruments are detailed in the table below:

LIKE Instruments	Contraindications
Drill	<p>These instruments are used for the installation of implants, so the contraindications of these instruments are the same as those of LIKE implants. It is contraindicated to place dental implants in patients with:</p> <ul style="list-style-type: none"> - medically unfit for an oral surgical procedure - with inadequate bone volume, unless an augmentation procedure can be considered - for whom it is impossible to place the adequate size or number of implants, or in the desirable position; treatment plan is then incompatible with the physiological functional constraints.
Implant driver (manual and contra angle)	
Drill extension	
Torque wrench	
Axis indicator	
Screwing key / Screwing	<p>These instruments are used for the installation of prostheses, the contraindications of these instruments are therefore the same as those of LIKE prosthetic surfaces.</p> <p>The prosthetic design must accommodate individual patient conditions such as bruxism or unfavourable jaw relationships to reduce the risk of overload or unusual wear, and treatment is contraindicated if adequate accommodation cannot be accomplished.</p>

INSTRUCTIONS FOR CLINICIAN

Each implant system has specific design characteristics for abutments and prosthetic components. Combining components that are not dimensioned or configured for correct mating can lead to mechanical failure of components, tissue damage or unsatisfactory aesthetic results.

We strongly recommend that clinicians, new as well as experimented implant users, always go through special training before undertaking a new treatment method. Always work with an experienced colleague the first time you employ a new treatment method.

EXPECTED PROFIT

The expected profit of the LIKE Implant System is the dental restoration

PATIENT INFORMATION:

The clinician is responsible for:

- informing the patient of the potential risks and adverse effects associated with the fitting of such an implant system,
- informing the patient that the effectiveness of the devices depends on his behaviour and his oral hygiene,
- ensuring that the patient reports him all suspicious movements between implant and abutment, and will have to undergo periodic inspections for detecting signs that precede abnormal functional symptoms.

PROCEDURAL PRECAUTIONS

The responsibility of the practitioner :

- To prevent the risk of inhalation and ingestion by the patient.
- To make sure to have the necessary ancillary for handling implant or prosthetic system.
- To check the screw torques, ancillaries functioning and retention force before use.
- To use the shortest drill as possible in accordance with the clinical requirements.
- To externally irrigate the surgical site while drilling.
- To not file contaminated instrument(s) with clean instrument(s) in the same container.

The drilling speed should not exceed 2000 rev / min.

The implant screwing speed should not exceed 25rev / min.

POTENTIAL ADVERSE EFFECTS

- Dislocation, rupture or loss of a component, infections, aesthetic problems, allergies to the different materials contained in prosthetic components.
- Aspirations of a component by the patient, allergy to any component
- Transient or permanent local nervous lesion, local sensitivity loss
- Bone fracture
- Infection

The surgeon is responsible for the complications that can result from an erroneous indication, of a faulty surgical technique, or a lack of sterility. In any case, these complications cannot be attributed to DEEP Company.

PACKAGING, STORAGE AND HANDLING

From delivery until use, store devices in a clean room away from products which could have a corrosive action.

In the event of an apparent packaging defect, the product must be returned to DEEP Company.

PROTOCOL OF CLEANING AND STERILIZATION BEFORE USE

Before first use, LIKE instruments must undergo an initial treatment. Indeed, the non-sterile equipment supplied by DEEP Company has been cleaned of manufacturing residues and then washed, but does not have a sufficient level of decontamination to permit it to be used directly. A decontamination and cleaning treatment is therefore essential prior to any sterilization.

The initial treatment protocol, validated by DEEP Company, is presented below. We recommend the use of detergents mentioned in protocol or an equivalent detergent, at neutral pH.

The choice of manual or automated cleaning is up to the surgeon.

1. Manual Cleaning or	<p>Equipment:</p> <ul style="list-style-type: none"> - Detergent (DENTASEPT ENZYMATIQUE d'Anios - Concentration 1%); - non-metallic soft brush; - tap water; - demineralised water (flushing) <p>Note: Be careful not to damage medical devices. Do not bang them together.</p> <ol style="list-style-type: none"> 1. Fully submerge the devices in a cleaning bath having a temperature of less than 50°C (Room temperature < 30°C). Any hinged devices should be fully opened and any multicomponent devices should be disassembled. Follow the soaking time recommended by the product manufacturer. Note: Unless specified otherwise, instruments and devices must be soaked for at least 15 minutes but avoid over soaking. 2. Brush the medical device inside and out, and any engraving, to remove persistent stains. Run the solution through device cavities to remove any residue. 3. Rinse thoroughly (Temp < 30°C) by soaking and/or spraying with water to remove detergent and soil; demineralised water should be used whenever possible.
1. Automated cleaning	<p>Equipment:</p> <ul style="list-style-type: none"> - Washer-disinfector according to NF EN ISO 15883-1&-2 (MIELE G 7825); - Detergent (Néodisher® MediClean - concentration 0.5%); - demineralised water. <ol style="list-style-type: none"> 1. Arrange the devices in a manner that allows water to flow into their cavities. 2. Start the cycle: <ul style="list-style-type: none"> - Prewash without product for 1 minute at room temperature - Wash: 40°C water + detergent - Bring temperature up to 50°C and maintain it for 5 minutes - Rinse with cold water for 1 minute at room temperature 3. Unload the instruments and inspect the surface and cavities of each device to make sure all visible soil has been removed. If needed, repeat the cycle or perform manual cleaning.
2. Automated disinfection	<p>If using a washer-disinfector, a final rinsing step with thermal disinfection (demineralised water temperature increased to 90°C and held for 5 minutes) can be added. The hold time depends on the size and power of the unit (ISO 15883-1).</p>
3. Drying	<p>Note: This step prevents microbial growth, removes any traces left from the washing and rinsing steps and optimises sterilisation effectiveness.</p> <ul style="list-style-type: none"> - If using the washer-disinfector, do not exceed 120°C (Drying for 20 min at 95°C). - Manual drying can be performed using a clean non-woven absorbent tissue or a clean lint-free cloth, and medical compressed air.

DEEP Company thanks you for your support and remains at your disposal for any additional information you may require

DEEP Company
ZA du bois Saint-Pierre
38280 JANNEYRIAS - FRANCE

Tel : +33 (0) 4 28 01 26 03
www.deep-company.com
contact@deep-company.com

IMPORTANT : Please read these instructions carefully before use

4. Packaging	Note: No matter which cleaning method is used, medical devices must be quickly placed into packages to avoid recontamination. Use materials suitable for sterilisation; the packaging used to ship the prosthesis systems cannot be sterilised. Make sure the pouch is large enough to hold the medical device, without stretching the edges. Do not place multiple devices in the same pouch.
5. Sterilisation	Autoclave with high pressure saturated steam at 134°C for at least 18 minutes (prion cycle). If several medical devices are being sterilised at once, make sure the autoclave is not overloaded.
6. Storage	Isolate the sterile packaged items so they do not come into contact with other devices.

MEANING OF THE SYMBOLS

The meaning of the symbols available on the external label is available in the following table :

SYMBOLS			
	Device manufacturer	LOT	Batch number of the Device
REF	Device reference		Have to not be used if the package has been damaged or opened
QTE	Indicates that the product is a medical device	CE 0123	CE Mark with the number of the Notify Body
MD	Indicates that the product is a medical device	UDI	Indicates a support that contains unique device identification
	Consult the instruction for use		Date of manufacture

DEVICE REPROCESSING

It is imperative to reprocess medical devices immediately after their use, in order to avoid the drying of dirt on reusable equipment. Prolonged contact with blood, pulp, bone, dentine, cement and other residues could damage the devices and compromise the effectiveness of the entire treatment.

1. Disassembly of the device, if necessary	Disassemble any medical devices having multiple parts.
2. Pre-disinfection	<ol style="list-style-type: none"> Remove any significant residue on medical devices by rinsing them with cold water or wiping them with lint-free tissue. Soak the device in an aldehyde-free detergent-disinfectant bath
3. Device reprocessing	Repeat steps 1 to 6 of the protocol from the section "Protocol of cleaning and sterilization before use"

Warning : Do not return soiled medical devices to their storage box

Reprocessing of the LIKE device should be performed according to the validated protocol below :

Warning : DEEP Company performed a reprocessing validation of up to 6 cycles. Thus, we strongly recommend that you do not exceed this number of cycle

The end-of-life for a reusable medical device is usually determined by the amount of wear and damage from a use or reprocessing different from the one recommend by DEEP Company.

SYSTEM COMPATIBILITY INFORMATION

The compatibility table for each range is given below. If in doubt, please contact us by email or telephone.

LIKE range	Compatible systems
LIKE-BS	<ul style="list-style-type: none"> - Implants and prosthetic parts of the LIKE-BS range; - Implants and prosthetic parts whose connection is similar to that of the LIKE-BS implant
LIKE-I	<ul style="list-style-type: none"> - Implants and prosthetic parts of the LIKE-I range; - Implants and prosthetic parts whose connection is similar to that of the LIKE-I implant
LIKE-R/CC	<ul style="list-style-type: none"> - Implants and prosthetic parts of the LIKE-CC range; - Prosthetic parts of the LIKE-R range - Implants and prosthetic parts whose connection is similar to that of the LIKE-CC implant
LIKE-A	<ul style="list-style-type: none"> - Implants and prosthetic parts of the LIKE-A range; - Implants and prosthetic parts whose connection is similar to that of the LIKE-A implant

PRODUCT RETURN

In the event of device return following an order error, it must be returned unused in its original unopened packaging.

All incidents related to the use of LIKE ancillaries must imperatively be reported to DEEP Company. Please clean and sterilize the device(s) according to the previous described procedures and send them in double sterile pouch. The return will be accompanied by the " Implant removal management" form available in the "documentary" section of our website.

ADDITIONNAL INFORMATION

For more information on the indication and use of different systems, please contact DEEP Company by the following email : contact@deep-company.com or on our website : www.deep-company.com

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