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**IMPORTANT : Please read these instructions carefully before use**

**LIABILITY**

The user of DEEP Company ancillary components has the duty to determine whether or not any product is suitable for the particular patient and circumstances. DEEP Company disclaims any liability, express or implied, and shall have no responsibility for any, direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of DEEP Company products. The clinician is also obliged to study the latest developments in regard to the DEEP Company products and their applications regularly. In cases of doubt, the clinician has to contact DEEP Company. Since the processing and surgical application of this product are under the control of the clinician, they are his/her responsibility. DEEP Company assumes no liability whatsoever for damages arising thereof.

**DESCRIPTION AND INTENDED USE**

DESIGN Ancillaries are medical devices allowing the placement of DESIGN Implants and prosthetics parts in the context of dental rehabilitations. The intended use of each instrument are given below :

Family product	Intended use
Drill	Drill into alveolar bone
Implant driver (manual and contra angle)	Take the implants in their packaging and place them in the alveolar bone
Drill extension	Facilitate the drill's use between theeth
Axis indicator	Control the axis of the implant and the depth of preparation
Screwing key / Screwing	Screw/ Unscrew the prosthetic components with the implant
Torque wrench	Control the tightening torque of the screws
Stop drills	To control precisely of the drilling depth during the preparation of the implant site

**Materials :**

Family product	Materials
Drill	Stainless steel
Implant driver (manual and contra angle)	
Drill extension	
Screwing key / Screwing	
Torque wrench	Titane biocompatible de grade 5
Axis indicator	

**INDICATIONS**

DESIGN instruments are used during dental surgery for dental rehabilitation in mandible and maxilla. These restorations range from a single replacement to the use of elements for overdenture applications as well as the full arch of bridgework.

**INTENDED PATIENT POPULATION**

DESIGN Implant System is intended for people with the absence of at least one tooth who does not have any contraindications identified in the instructions for use.

For children, the placement of implants is not recommended until the end of jaw growth has been duly established. Preoperative gum and bone deficiencies can compromise the aesthetic result or result in unfavourable implant positioning.

**CONTRAINDICATIONS**

The contraindications associated with DESIGN instruments are the same as those indicated in the instruction for use of the DESIGN implants and prosthetics components.

**INSTRUCTIONS FOR CLINICIAN**

Each implant system has specific design characteristics for abutments and prosthetic components. Combining components that are not dimensioned or configured for correct mating can lead to mechanical failure of components, tissue damage or unsatisfactory aesthetic results.

We strongly recommend that clinicians, new as well as experimented implant users, always go through special training before undertaking a new treatment method. Always work with an experienced colleague the first time you employ a new treatment method.

**EXPECTED PROFIT**

The expected profit of the DESIGN Implant System is the dental restoration

**PATIENT INFORMATION :**

The clinician is responsible for:

- informing the patient of the potential risks and adverse effects associated with the fitting of such an implant system,
- informing the patient that the effectiveness of the devices depends on his behaviour and his oral hygiene,
- ensuring that the patient reports him all suspicious movements between implant and abutment, and will have to undergo periodic inspections for detecting signs that precede abnormal functional symptoms.

**PROCEDURAL PRECAUTIONS**

The responsibility of the practitioner :

- To prevent the risk of inhalation and ingestion by the patient.
  - To make sure to have the necessary ancillary for handling implant or prosthetic system.
  - To check the screw torques, ancillaries functioning and retention force before use.
  - To use the shortest drill as possible in accordance with the clinical requirements.
  - To externally irrigate the surgical site while drilling.
  - To not file contaminated instrument(s) with clean instrument(s) in the same container.
- The drilling speed should not exceed 2000 rev / min.  
The implant screwing speed should not exceed 25rev / min.

**POTENTIAL ADVERSE EFFECTS**

- Dislocation, rupture or loss of a component, infections, aesthetic problems, allergies to the different materials contained in prosthetic components.
  - Aspirations of a component by the patient, allergy to any component
  - Transient or permanent local nervous lesion, local sensitivity loss
  - Bone fracture
  - Infection
- The surgeon is responsible for the complications that can result from an erroneous indication, of a faulty surgical technique, or a lack of sterility. In any case, these complications cannot be attributed to DEEP Company.

**PACKAGING, STORAGE AND HANDLING**

From delivery until use, store devices in a clean room away from products which could have a corrosive action.

In the event of an apparent packaging defect, the product must be returned to DEEP Company.

**PROTOCOL OF CLEANING AND STERILIZATION BEFORE FIRST USE**

The non-sterile instruments supplied by DEEP Company have been cleaned of manufacturing residues and then washed but do not have a sufficient level of decontamination to allow them to be used directly. Therefore, before first use, all DESIGN instruments may be cleaned, disinfected and sterilized according to the protocol mentioned below.

The protocol was drawn up in accordance with the SF2S guide to Good Sterilisation Practice for Reusable Medical Devices, updated in 2021.

<b>Concerned devices</b>	All DESIGN ancillaries are covered by this manual
<b>Warning</b>	The cleaning, disinfection and sterilisation protocols have been validated by DEEP Company. We therefore recommend that you follow each step scrupulously (parameters, detergent, etc.). For cleaning and disinfection, DEEP Company has used the detergent Neodisher® MediClean Dental (DrWeiger) for the validation of the protocols. However, other products may be used if they are CE marked and aldehyde-free. The use of other detergents is the responsibility of the user and must be validated before use.  Protective clothing (mask, gloves and goggles) must be worn when handling cleaning and disinfecting detergent.

Instructions																			
<b>1. Preparation before cleaning</b>	<p>1. Remove the instruments from the surgical kit and wash them separately;</p> <p>2. Disassemble the following instruments :</p> <ul style="list-style-type: none"> <li>- for rotary instruments, separate the tool (drill) from the hand tool in accordance with the latter's instructions for use;</li> <li>- for manual screwdrivers, please separate them from their screw wrenches (torque or ratchet) according to the instructions for the latter;</li> <li>- for multi abutment manual screwdriver (ref : PDF30-i94040001), please follow the dismantling instructions below:               <ol style="list-style-type: none"> <li>1. Pull the knob upwards, then unscrew it;</li> <li>2. Remove the knob from the rest of the spanner body so that it can be cleaned.</li> </ol> </li> </ul>																		
	 <p>With the exception of the instruments listed above, the other instruments cannot be disassembled.</p>																		
<b>2. Cleaning &amp; disinfection : Automated</b>	<p><b>Equipment :</b> For this operation, DEEP Company recommends automatic cleaning and disinfection of its dental instruments using a washer-disinfector that complies with standard NF EN ISO 15883-1&amp;-2. The washer/disinfector must be validated beforehand by the user and regularly maintained.</p> <p><b>Detergent :</b> DEEP Company used Neodisher® MediClean Dental detergent (DrWeiger) at a concentration of 0.2%. The concentration used can vary from 0.2% to 1% depending on the degree of soiling.</p> <p><b>Steps :</b></p> <ol style="list-style-type: none"> <li>1. Place disassembled devices in a wash basket suitable for the washer-disinfector and the configuration of the instruments. Orient instruments in manner that allows water to flow into their cavities. Ensure that all surfaces are exposed to the cleaning solution.</li> <li>⚠ Do not overload the wash basket. The load is at the user's initiative and must be validated by the user during the initial qualification of the equipment.</li> </ol> <p>2. Start the next cycle :</p>																		
	<table border="1"> <thead> <tr> <th>Steps</th> <th>Water*</th> <th>Minimal temperature</th> <th>Minimal time</th> </tr> </thead> <tbody> <tr> <td><b>Prewash</b></td> <td>Tap water</td> <td>Cold water : 15°C – 30°C</td> <td>1 minute</td> </tr> <tr> <td><b>Wahsing</b></td> <td rowspan="2">Demineralised water + detergent</td> <td>Heated water : 40-60°C</td> <td>10 minutes</td> </tr> <tr> <td><b>Thermal disinfection</b></td> <td>Heated water : 85°C-90°C</td> <td>5 minutes</td> </tr> <tr> <td><b>Rinsing</b></td> <td>Demineralised water</td> <td>Cold water : 15°C – 30°C</td> <td>1 minute</td> </tr> </tbody> </table>	Steps	Water*	Minimal temperature	Minimal time	<b>Prewash</b>	Tap water	Cold water : 15°C – 30°C	1 minute	<b>Wahsing</b>	Demineralised water + detergent	Heated water : 40-60°C	10 minutes	<b>Thermal disinfection</b>	Heated water : 85°C-90°C	5 minutes	<b>Rinsing</b>	Demineralised water	Cold water : 15°C – 30°C
Steps	Water*	Minimal temperature	Minimal time																
<b>Prewash</b>	Tap water	Cold water : 15°C – 30°C	1 minute																
<b>Wahsing</b>	Demineralised water + detergent	Heated water : 40-60°C	10 minutes																
<b>Thermal disinfection</b>		Heated water : 85°C-90°C	5 minutes																
<b>Rinsing</b>	Demineralised water	Cold water : 15°C – 30°C	1 minute																

**DEEP Company thanks you for your support and remains at your disposal for any additional information you may require**

DEEP Company  
ZA du bois Saint-Pierre  
38280 JANNEYRIAS - FRANCE

Tel : +33 (0) 4 28 01 26 03  
[www.deep-company.com](http://www.deep-company.com)  
contact@deep-company.com

### IMPORTANT : Please read these instructions carefully before use

	Drying	N/A	95°C	20 minutes									
	<p>* Water quality specifications for cleaning and rinsing are available at the end of the protocol.</p> <p>3. When unloading, inspect the surfaces and cavities of medical devices to verify complete removal of visible dirt. If necessary, repeat the cycle or carry out manual cleaning.</p> <p>4. Likewise, check that the instruments do not show any deformation, damage or defects (loss of marking), caused by the cleaning step, and which could affect the safety and performance of the device. If this is the case, please stop the operation and return the defective products to DEEP Company (see "material return" section).</p>												
3. Drying	<p>The drying step is integrated into the washing/disinfection step above. Visually check that the devices are dry when unloading the devices. If this is not the case, carefully wipe all device surfaces with a clean, lint-free cloth (single use) until the instruments are dry.</p>												
4. Packaging	<p><b>Type of packaging to be used :</b> DEEP Company recommends the use of a pre-formed sterile barrier system (sterilisation sheath) that complies with standard NF EN ISO 11 607-1 and is suitable for autoclave sterilisation. In addition, these sheaths must be adapted to the dimensions of the surgical kit.</p> <p>For validation of sterilization, DEEP Company used packaging consisting of two sheaths: an internal sheath (reference: SCL10152) serving as a preformed sterile barrier, and external sheath (reference: SCLI 2182) to protect the sterile barrier, all two from the brand "Crosstex".</p> <p>▲ Please check the condition of the sterilisation sheath before packaging.</p> <p><b>Steps :</b></p> <ol style="list-style-type: none"> <li>For multi abutment manual screwdriver (ref : PDF30-i94040001), please assemble in reverse order to disassembly instructions ;</li> <li>Place the instruments in the appropriate surgical kit supplied by DEEP Company.</li> <li>To avoid recontamination, pack the surgical kit in the sterilisation sheath in accordance with its instructions for use.</li> <li>Check that the packaging is properly closed, then affix a traceability label containing the following information: the product reference and IUD, the type of sterilisation to be used, the use-by date (depending on the sachet used) and the applicable storage conditions if applicable.</li> </ol>												
5. Sterilization	<p><b>Equipment :</b> DEEP Company recommends performing the sterilisation of its dental instruments using a CE-marked, validated autoclave that complies with standard NF EN ISO 17 665. The autoclave must be regularly maintained in order to maintain sterilization efficiency.</p> <p><b>Steps :</b></p> <ol style="list-style-type: none"> <li>Place the surgical kit, previously packaged, in the autoclave basket to ensure even steam penetration through the packaging. Make sure there is no contact between the packaging and the wall of the steriliser ;</li> </ol> <p>▲ Do not overload the autoclave basket. The load is the user's responsibility and must be validated by the user during the initial qualification of the equipment.</p> <ol style="list-style-type: none"> <li>To check the effectiveness of the sterilisation cycle, place a physical-chemical indicator on the load to be sterilised.</li> <li>Start a sterilisation cycle using the parameters below * :</li> </ol> <table border="1"> <thead> <tr> <th>Countries</th> <th>Sterilisation parameters</th> <th>Drying time</th> </tr> </thead> <tbody> <tr> <td>France</td> <td>134°C (- 0°C / + 3°C) pendant 18 minutes</td> <td>30 min</td> </tr> <tr> <td>Europe</td> <td>134°C (- 0°C / + 3°C) pendant 3 minutes</td> <td>30 min</td> </tr> </tbody> </table> <p>* the choice of cycle must comply with the standards and regulations in force in the country of sale</p> <ol style="list-style-type: none"> <li>After unloading the load, check the physical-chemical indicator (colour change of sterilisation indicators) and the sterilisation efficiency (integrity of the packaging, no humidity, integrity of the traceability label, digital recordings of the cycle parameters). Keep the indicator and the checks in the patient file to ensure traceability.</li> </ol>				Countries	Sterilisation parameters	Drying time	France	134°C (- 0°C / + 3°C) pendant 18 minutes	30 min	Europe	134°C (- 0°C / + 3°C) pendant 3 minutes	30 min
Countries	Sterilisation parameters	Drying time											
France	134°C (- 0°C / + 3°C) pendant 18 minutes	30 min											
Europe	134°C (- 0°C / + 3°C) pendant 3 minutes	30 min											
6. Storage	<p>Please keep the instruments in their sterilization packaging in a dry and clean environment until the expiration date (inscribed on the traceability label). Please consult the notice of packaging used for further information on storage conditions, if applicable.</p> <p>Before use, please check the packaging and instruments (integrity of the packaging, sign of humidity, etc.). In the event of damage to the packaging, sterility cannot be guaranteed.</p>												

Proprieties	Limit values
Conductivity	< 1.3µS/cm à 25°C
Bacteria	< 10 CFU/100 mL
Endotoxin	<0.25 EU/mL

#### DEVICE REPROCESSING

This protocol is applicable after the first use of the instruments and after each use.

<b>Concerned devices</b>	All DESIGN instruments are covered by these instructions for use
<b>Warning</b>	This operation involves handling devices contaminated with blood, bones and saliva. It is therefore essential that the person in charge of this operation wears appropriate protective clothing (medical mask, gloves, goggles and single-use impermeable apron) to prevent microbiological and chemical risks.
<b>Limits to reprocessing</b>	The lifetime of DESIGN instruments has been validated at 6 uses.  However, if the instruments show any deformation, damage or defects (loss of marking) that may affect the safety and performance of the device, please return the defective products to DEEP Company.

DESIGN devices must be reprocessed in accordance with the protocol validated below. The protocol has been drawn up in accordance with the SF2S guide to Good Sterilisation Practice for Reusable Medical Devices, updated in 2021.

<b>1. Disassembly of devices</b>	For devices that can be disassembled, disassembly them as indicated in section "1. Preparation before cleaning" in the previous paragraph.  <b>Equipment :</b> DEEP Company recommends : - the use of a soft brush with nylon bristles, adapted to the size of the instruments. The use of wire brushes or scouring pads is prohibited to avoid scratching the surface of the instruments; - a soaking tank made of corrosion-resistant material and adapted to the size of the instruments..  <b>Detergent :</b> DEEP Company used Neodisher® MediClean Dental detergent (DrWeigert) at a concentration of 0.2%. The concentration used can vary from 0.2% to 1% depending on the degree of soiling.  <b>Steps :</b>  1. Soak the instruments in a soaking tank containing a detergent solution. The water temperature should be between 25°C and 40°C. The soaking time should not be less than 10 minutes and should not exceed 30 minutes to avoid corrosion of the instruments.;
<b>2. Initial treatment at the point of use</b>	2. Remove large residues from the instruments using a soft brush until the instruments are visually free of all traces of dirt (blood, paste, cement, etc.);  3. Rinse dental instruments with demineralised water at a temperature of between 15°C and 30°C for at least one minute to remove the detergent;  4. Carefully wipe all device surfaces with a clean lint-free cloth (single use) until the instruments are dry.  5. Pre-cleaned devices should be placed in suitable, lidded medical transport container to reduce the risk of exposure of employees and the environment to pathogens. The container can be transported by hand or using a suitable transport trolley.  6. To avoid any risk of recontamination, DEEP Company recommends that pre-cleaned devices should not be stored, and that mechanical cleaning of instruments should be carried out within a maximum of 30 minutes.
<b>3. Reprocessing of devices</b>	Repeat steps 2 to 6 of the protocol in the "Cleaning & sterilisation protocol before first use" section.

#### SYSTEM COMPATIBILITY INFORMATION

DESIGN instruments are compatible with DESIGN implants and prosthetic components.

#### PRODUCT RETURN

In the event of device return following an order error, it must be returned unused in its original unopened packaging.

All incidents related to the use of DESIGN ancillaries must imperatively be reported to DEEP Company. Please clean and sterilize the device(s) according to the previous described procedures and send them in double sterile pouch. The return will be accompanied by the " Implant removal management" form available in the "documentary" section of our website.

#### ADDITIONNAL INFORMATION

For more information on the indication and use of different systems, please contact DEEP Company by the following email : [contact@deep-company.com](mailto:contact@deep-company.com) or on our website : [www.deep-company.com](http://www.deep-company.com)


#### MEANING OF THE SYMBOLS

The meaning of the symbols available on the external label is available in the following table :

**Note :** The practitioner should carry out checks on the water used for manual and automatic cleaning. DEEP Company recommends using tap water or demineralised water that meets the following local pharmacopoeia specifications:

**IMPORTANT : Please read these instructions carefully before use**

<b>SYMBOLS</b>		Device manufacturer		Batch number of the Device		Device reference
		Single use only		CE Mark with the number of the Notify Body		Indicates that the product is a medical device
		Consult the instruction for use		Date of manufacture		Indicates a support that contains unique device identification
		Have to not be used if the package has been damaged or opened		Indicates that the product is a medical device		